

N4BP2L2 Antibody



PACO02316

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

N4BP2L2

Uniprot

Q92802

Synonyms:

N4BP2L2; CG005; PFAAP5; NEDD4-binding protein 2-like 2; Phosphonoformate immuno-associated protein 5


Immunogen:

Synthesized peptide derived from the Internal region of human PFAAP5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of LOVO cells using PFAAP5 Polyclonal Antibody.