

RGS14 Antibody



PACO02349

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

RGS14

Uniprot

O43566

Synonyms:

RGS14; Regulator of G-protein signaling 14; RGS14


Immunogen:

Synthesized peptide derived from the Internal region of human RGS14.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 3T3 cells using RGS14 Polyclonal Antibody.