

RRH Antibody



PACO02355

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

RRH

Uniprot

O14718

Synonyms:

RRH; Visual pigment-like receptor peropsin


Immunogen:

Synthesized peptide derived from the Internal region of human Peropsin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Peropsin Polyclonal Antibody.