

# CYSLTR1 Antibody



PACO02368

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Monkey

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:****Gene ID:**

CYSLTR1

**Uniprot**

Q9Y271

**Synonyms:**

CYSLTR1; CYSLT1; Cysteinyl leukotriene receptor 1; CysLTR1; Cysteinyl leukotriene D4 receptor; LTD4 receptor; G-protein coupled receptor HG55; HMTMF81

**Immunogen:**

Synthesized peptide derived from the Internal region of human CysLTR1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of COS-7 cells using CysLTR1 Polyclonal Antibody.

