

PNPLA6 Antibody



PACO02379

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

PNPLA6

Uniprot

Q8IY17

Synonyms:

PNPLA6; NTE; Neuropathy target esterase; Patatin-like phospholipase domain-containing protein 6


Immunogen:

Synthesized peptide derived from the N-terminal region of human NTE.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of SH-SY5Y cells using NTE Polyclonal Antibody.