

ATXN2L Antibody



PACO02416

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ATXN2L

Uniprot

Q8WWM7

Synonyms:

ATXN2L; A2D; A2LG; A2LP; A2RP; Ataxin-2-like protein; Ataxin-2 domain protein; Ataxin-2-related protein

Immunogen:

Synthesized peptide derived from the Internal region of human Ataxin-2L.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293T cells using Ataxin-2L Polyclonal Antibody.

