

Phospho-TIRAP (Y86) Antibody



PACO02425

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

TIRAP

Uniprot

P58753

Synonyms:

TIRAP; MAL; Toll/interleukin-1 receptor domain-containing adapter protein; TIR domain-containing adapter protein; Adaptor protein Wyatt; MyD88 adapter-like protein


Immunogen:

Synthesized peptide derived from human TIRAP around the phosphorylation site of Y86.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of K562 using p-TIRAP (Y86) antibody.