

CLN5 Antibody



PACO02444

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

CLN5

Uniprot

O75503

Synonyms:

CLN5; Ceroid-lipofuscinosis neuronal protein 5; Protein CLN5


Immunogen:

Synthesized peptide derived from the Internal region of human CLN5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293T cells using CLN5 Polyclonal Antibody.