

CLPTM1 Antibody



PACO02447

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

CLPTM1

Uniprot

O96005

Synonyms:

CLPTM1; Cleft lip and palate transmembrane protein 1


Immunogen:

Synthesized peptide derived from the Internal region of human CLPTM1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS-7 cells using CLPTM1 Polyclonal Antibody.