

FOXN4 Antibody



PACO02453

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

FOXN4

Uniprot

Q96NZ1

Synonyms:

FOXN4; Forkhead box protein N4


Immunogen:

Synthesized peptide derived from the C-terminal region of human FoxN4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of LOVO cells using FoxN4 Polyclonal Antibody.