

OR10Z1 Antibody



PACO02474

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR10Z1

Uniprot

Q8NGY1

Synonyms:

OR10Z1; Olfactory receptor 10Z1; Olfactory receptor OR1-15

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 10Z1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of COS-7 cells using Olfactory receptor 10Z1 Polyclonal Antibody.

