

Phospho-CSF1R (Y809) Antibody



PACO02483

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

CSF1R

Uniprot

P07333

Synonyms:

CSF1R; FMS; Macrophage colony-stimulating factor 1 receptor; CSF-1 receptor; CSF-1-R; CSF-1R; M-CSF-R; Proto-oncogene c-Fms; CD antigen CD115

Immunogen:

Synthesized peptide derived from human c-Fms around the phosphorylation site of Y809.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western blot analysis of 293 using p-c-Fms (Y809) antibody.

