

Phospho-CSK (S364) Antibody



PACO02489

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

CSK

Uniprot

P41240

Synonyms:

CSK; Tyrosine-protein kinase CSK; C-Src kinase; Protein-tyrosine kinase CYL


Immunogen:

Synthesized peptide derived from human Csk around the phosphorylation site of S364.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of VEC cells using Phospho-Csk (S364) Polyclonal Antibody.