

# Phospho-CSNK1A1 (Y321) Antibody



PACO02491

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

CSNK1A1

**Uniprot**

P48729

**Synonyms:**

CSNK1A1; Casein kinase I isoform alpha; CKI-alpha; CK1

**Immunogen:**

Synthesized peptide derived from human Casein Kinase I alpha around the phosphorylation site of Y321.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of VEC cells using Phospho-Casein Kinase Ialpha (Y321) Polyclonal Antibody.

