

CMTM2 Antibody



PACO02501

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CMTM2

Uniprot

Q8TAZ6

Synonyms:

CMTM2; CKLFSF2; CKLF-like MARVEL transmembrane domain-containing protein 2; Chemokine-like factor superfamily member 2

Immunogen:

Synthesized peptide derived from the C-terminal region of human CMTM2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of K562 cells using CMTM2 Polyclonal Antibody.

