

# CST1 Antibody



PACO02504

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

CST1

**Uniprot**

P01037

**Synonyms:**

CST1; Cystatin-SN; Cystain-SA-I; Cystatin-1; Salivary cystatin-SA-1

**Immunogen:**

Synthesized peptide derived from the Internal region of human Cystatin SN.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of Jurkat cells using Cystatin SN Polyclonal Antibody.