

ZNF420 Antibody



PACO02509

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ZNF420

Uniprot

Q8TAQ5

Synonyms:

ZNF420; Zinc finger protein 420

Immunogen:

Synthesized peptide derived from the Internal region of human ZNF420.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of K562 cells using ZNF420 Polyclonal Antibody.

