

CTBP1 Antibody



PACO02513

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

CTBP1

Uniprot

Q13363

Synonyms:

CTBP1; CTBP; C-terminal-binding protein 1; CtBP1

Immunogen:

Synthesized peptide derived from human CtBP1 around the non-phosphorylation site of S422.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of Jurkat cells using CtBP1 Polyclonal Antibody.

