

Phospho-CTNNB1 (Y654) Antibody



PACO02517

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

CTNNB1

Uniprot

P35222

Synonyms:

CTNNB1; CTNNB; OK/SW-cl.35; Catenin beta-1; Beta-catenin


Immunogen:

Synthesized peptide derived from human Catenin-Beta around the phosphorylation site of Y654.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using Phospho-Catenin-Beta (Y654) Polyclonal Antibody.