

# CTNNB1 Antibody



PACO02518

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

CTNNB1

**Uniprot**

P35222

**Synonyms:**

CTNNB1; CTNNB; OK/SW-cl.35; Catenin beta-1; Beta-catenin

**Immunogen:**

Synthesized peptide derived from human Catenin-Beta around the non-phosphorylation site of S552.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of A549 cells using Catenin-Beta Polyclonal Antibody.