

CTSD Antibody



PACO02523

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

CTSD

Uniprot

P07339

Synonyms:

CTSD; CPSD; Cathepsin D


Immunogen:

Synthesized peptide derived from the Internal region of human Cathepsin D.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Cathepsin D Polyclonal Antibody.