

GDF7 Antibody



PACO02526

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

GDF7

Uniprot

Q7Z4P5

Synonyms:

GDF7; Growth/differentiation factor 7; GDF-7

Immunogen:

Synthesized peptide derived from the Internal region of human GDF-7.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HT-29 cells using GDF-7 Polyclonal Antibody.

