

SLC16A14 Antibody



PACO02527

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

SLC16A14

Uniprot

Q7RTX9

Synonyms:

SLC16A14; MCT14; Monocarboxylate transporter 14; MCT 14; Solute carrier family 16 member 14

Immunogen:

Synthesized peptide derived from the Internal region of human MCT14.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HEPG2 cells using MCT14 Polyclonal Antibody.

