

PACO02532

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

ARSK

Uniprot

Q6UWY0

Synonyms:

ARSK; TSULF; Arylsulfatase K; ASK; Telethon sulfatase


Immunogen:

Synthesized peptide derived from the C-terminal region of human Arylsulfatase K.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using Arylsulfatase K Polyclonal Antibody.