

CYP2C8 Antibody



PACO02537

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CYP2C8

Uniprot

P10632

Synonyms:

CYP2C8; Cytochrome P450 2C8; CYP11C8; Cytochrome P450 11C2; Cytochrome P450 MP-12; Cytochrome P450 MP-20; Cytochrome P450 form 1; S-mephenytoin 4-hydroxylase


Immunogen:

Synthesized peptide derived from the Internal region of human CYP2C8.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using CYP2C8 Polyclonal Antibody.