

CYP2C8/CYP2C9/CYP2C18/CYP2C19 Antibody



PACO02540

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

CYP2C8/CYP2C9/CYP2C18/CYP2C19

Uniprot

P10632/P11712/P33260/P33261

Synonyms:

CYP2C8; Cytochrome P450 2C8; CYP11C8; Cytochrome P450 IIC2; Cytochrome P450 MP-12; Cytochrome P450 MP-20; Cytochrome P450 form 1; S-mephenytoin 4-hydroxylase; CYP2C9; CYP2C10; Cytochrome P450 2C9; (R)-limonene 6-monooxygenase; (S)-limonene


Immunogen:

Synthesized peptide derived from the Internal region of human CYP2C8/9/18/19.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of LOVO cells using CYP2C8/9/18/19 Polyclonal Antibody.