

Phospho-SGK223 (Y413) Antibody



PACO02543

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

PRAG1

Uniprot

Q86YV5

Synonyms:

SGK223; Tyrosine-protein kinase Sgk223; Sugen kinase 223

Immunogen:

Synthesized peptide derived from human Pragmin around the phosphorylation site of Y413.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of A549 cells using Phospho-Pragmin (Y413) Polyclonal Antibody.