

NKX6-3 Antibody



PACO02546

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

NKX6-3

Uniprot

A6NJ46

Synonyms:

NKX6-3; Homeobox protein Nkx-6.3

Immunogen:

Synthesized peptide derived from the Internal region of human Nkx-6.3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of K562 cells using Nkx-6.3 Polyclonal Antibody.

