

CYP17A1 Antibody



PACO02549

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

CYP17A1

Uniprot

P05093

Synonyms:

CYP17A1; CYP17; S17AH; Steroid 17-alpha-hydroxylase/17; 20 lyase; CYPXVII;
Cytochrome P450 17A1; Cytochrome P450-C17; Cytochrome P450c17; Steroid 17-
alpha-monooxygenase

Immunogen:

Synthesized peptide derived from the Internal region of human CYP17A1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293 cells using CYP17A1 Polyclonal Antibody.

