

ZDHHC15 Antibody



PACO02551

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

ZDHHC15

Uniprot

Q96MV8

Synonyms:

ZDHHC15; Palmitoyltransferase ZDHHC15; Zinc finger DHHC domain-containing protein 15; DHHC-15


Immunogen:

Synthesized peptide derived from the C-terminal region of human DHHC-15.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using DHHC-15 Polyclonal Antibody.