

# GPR180 Antibody



PACO02557

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:****Gene ID:**

GPR180

**Uniprot**

Q86V85

**Synonyms:**

GPR180; ITR; Integral membrane protein GPR180; Intimal thickness-related receptor

**Immunogen:**

Synthesized peptide derived from the Internal region of human GPR180.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of LOVO cells using GPR180 Polyclonal Antibody.

