

Phospho-DAXX (S668) Antibody



PACO02563

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

DAXX

Uniprot

Q9UER7

Synonyms:

DAXX; BING2; DAP6; Death domain-associated protein 6; Daxx; hDaxx; ETS1-associated protein 1; EAP1; Fas death domain-associated protein

Immunogen:

Synthesized peptide derived from human Daxx around the phosphorylation site of S668.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293 cells using Phospho-Daxx (S668) Polyclonal Antibody.

