

ADAMTS19 Antibody



PACO02584

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

ADAMTS19

Uniprot

Q8TE59

Synonyms:

ADAMTS19; A disintegrin and metalloproteinase with thrombospondin motifs 19; ADAM-TS 19; ADAM-TS19; ADAMTS-19


Immunogen:

Synthesized peptide derived from the Internal region of human ADAMTS-19.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using ADAMTS-19 Polyclonal Antibody.