

Phospho-DNM1 (S778) Antibody



PACO02590

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

DNM1

Uniprot

Q05193

Synonyms:

DNM1; DNM; Dynamin-1


Immunogen:

Synthesized peptide derived from human Dynamin I around the phosphorylation site of S778.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Mouse-Brain cells using Phospho-Dynamin I (S778) Polyclonal Antibody.