

DOK1 Antibody



PACO02594

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

DOK1

Uniprot

Q99704

Synonyms:

DOK1; Docking protein 1; Downstream of tyrosine kinase 1; p62(dok); pp62

Immunogen:

Synthesized peptide derived from human Dok-1 around the non-phosphorylation site of Y362.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of Jurkat cells using Dok-1 Polyclonal Antibody.

