

Phospho-EEF2 (T56) Antibody



PACO02614

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

EEF2

Uniprot

P13639

Synonyms:

EEF2; EF2; Elongation factor 2; EF-2

Immunogen:

Synthesized peptide derived from human EF-2 around the phosphorylation site of T56.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of SKOV3 cells using Phospho-EF-2 (T56) Polyclonal Antibody.

