

Phospho-EFNB1 (Y317) Antibody



PACO02615

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

EFNB1

Uniprot

P98172

Synonyms:

EFNB1; EFL3; EPLG2; LERK2; Ephrin-B1; EFL-3; ELK ligand; ELK-L; EPH-related receptor tyrosine kinase ligand 2; LERK-2

Immunogen:

Synthesized peptide derived from human Ephrin-B1 around the phosphorylation site of Y317.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of mouse cells using Phospho-Ephrin-B1 (Y317) Polyclonal Antibody.

