

Phospho-EPHA2/EPHA5 (Y594) Antibody



PACO02631

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

EPHA2/EPHA5

Uniprot

P29317/P54756

Synonyms:

EPHA2; ECK; Ephrin type-A receptor 2; Epithelial cell kinase; Tyrosine-protein kinase receptor ECK; EPHA5; BSK; EHK1; HEK7; TYRO4; Ephrin type-A receptor 5; Brain-specific kinase; EPH homology kinase 1; EHK-1; EPH-like kinase 7; EK7; hEK7


Immunogen:

Synthesized peptide derived from human EphA2/5 around the phosphorylation site of Y594.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of JK cells using Phospho-EphA2/5 (Y594) Polyclonal Antibody.