

APOBEC3A Antibody



PACO02638

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

APOBEC3A

Uniprot

P31941

Synonyms:

APOBEC3A; Probable DNA dC->dU-editing enzyme APOBEC-3A; Phorbolin-1

Immunogen:

Synthesized peptide derived from the Internal region of human APOBEC3A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of Jurkat cells using APOBEC3A Polyclonal Antibody.

