

FOXD4L1 Antibody



PACO02639

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

FOXD4L1

Uniprot

Q9NU39

Synonyms:

FOXD4L1; Forkhead box protein D4-like 1; FOXD4-like 1


Immunogen:

Synthesized peptide derived from the C-terminal region of human FoxD4L1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using FoxD4L1 Polyclonal Antibody.