

SLC16A13 Antibody



PACO02644

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

SLC16A13

Uniprot

Q7RTY0

Synonyms:

SLC16A13; MCT13; Monocarboxylate transporter 13; MCT 13; Solute carrier family 16 member 13


Immunogen:

Synthesized peptide derived from the C-terminal region of human MCT13.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 3T3 cells using MCT13 Polyclonal Antibody.