

OR5A1 Antibody



PACO02700

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

OR5A1

Uniprot

Q8NGJ0

Synonyms:

OR5A1; OR5A1P; Olfactory receptor 5A1; OST181; Olfactory receptor OR11-249


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 5A1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MCF-7 cells using Olfactory receptor 5A1 Polyclonal Antibody.