

OR4D6 Antibody



PACO02701

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

OR4D6

Uniprot

Q8NGJ1

Synonyms:

OR4D6; Olfactory receptor 4D6; Olfactory receptor OR11-250


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4D6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of LOVO cells using Olfactory receptor 4D6 Polyclonal Antibody.