

CEP152 Antibody



PACO02742

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

CEP152

Uniprot

O94986

Synonyms:

CEP152; KIAA0912; Centrosomal protein of 152 kDa; Cep152


Immunogen:

Synthesized peptide derived from the Internal region of human CEP152.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using CEP152 Polyclonal Antibody.