

GATA3/GATA2 Antibody



PACO02777

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

GATA3/GATA2

Uniprot

P23771/P23769

Synonyms:

GATA3; Trans-acting T-cell-specific transcription factor GATA-3; GATA-binding factor 3; GATA2; Endothelial transcription factor GATA-2; GATA-binding protein 2


Immunogen:

Synthesized peptide derived from human GATA-3 around the non-phosphorylation site of S308.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of Mouse-brain lysis using GATA-2/3 antibody.