

Phospho-PDCD4 (S457) Antibody



PACO02785

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

PDCD4

Uniprot

Q53EL6

Synonyms:

PDCD4; H731; Programmed cell death protein 4; Neoplastic transformation inhibitor protein; Nuclear antigen H731-like; Protein 197/15a


Immunogen:

Synthesized peptide derived from human Pdc4 around the phosphorylation site of S457.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 cells using Phospho-Pdcd-4 (S457) Polyclonal Antibody.