

Phospho-GRIN2B (Y1336) Antibody



PACO02797

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

GRIN2B

Uniprot

Q13224

Synonyms:

GRIN2B; NMDAR2B; Glutamate [NMDA] receptor subunit epsilon-2; N-methyl D-aspartate receptor subtype 2B; NMDAR2B; NR2B; N-methyl-D-aspartate receptor subunit 3; NR3; hNR3


Immunogen:

Synthesized peptide derived from human NMDAepsilon2 around the phosphorylation site of Y1336.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293T-UV cells using Phospho-NMDAepsilon2 (Y1336) Polyclonal Antibody.