

Phospho-ARHGAP35 (Y1087) Antibody



PACO02800

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

ARHGAP35

Uniprot

Q9NRY4

Synonyms:

ARHGAP35; GRF1; GRLF1; KIAA1722; Rho GTPase-activating protein 35; Glucocorticoid receptor DNA-binding factor 1; Glucocorticoid receptor repression factor 1; GRF-1; Rho GAP p190A; p190-A


Immunogen:

Synthesized peptide derived from human GRF-1 around the phosphorylation site of Y1087.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293T cells using Phospho-GRF-1 (Y1087) Polyclonal Antibody.