

Phospho-Histone H1 (T17) Antibody



PACO02807

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

HIST1H1B/HIST1H1D/HIST1H1E

Uniprot

P16401/P16402/P10412

Synonyms:

HIST1H1B; H1F5; Histone H1.5; Histone H1a; Histone H1b; Histone H1s-3; HIST1H1D;
H1F3; Histone H1.3; Histone H1c; Histone H1s-2; HIST1H1E; H1F4; Histone H1.4;
Histone H1b; Histone H1s-4

Immunogen:

Synthesized peptide derived from human Histone H1 around the phosphorylation site of T17.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western blot analysis of SKOV3 using p-Histone H1 (T17) antibody.

