

Phospho-ACACA (S80) Antibody



PACO02811

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

ACACA

Uniprot

Q13085

Synonyms:

ACACA; ACAC; ACC1; ACCA; Acetyl-CoA carboxylase 1; ACC1; ACC-alpha

Immunogen:

Synthesized peptide derived from human ACCalpha around the phosphorylation site of S80.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western blot analysis of 3T3 using p-ACCalpha (S80) antibody.

